



## Instructions Manual

### Full HD 1080p SD-Card Image Capture Device

### Firmware Update v1.30

This firmware update is for products OP-009 080 and OP-009 081. Besides some bug fixes, most important improvements in the updated firmware v1.30 are:

- Image capture speed improved – time to save image 3.5 - 4 sec\* (it was ~ 7 sec\*)
- It is possible to save image with graticule in 7 - 8 sec\* (OP-009 081)

#### Installation:

1. Copy the file Optilia\_firmware.oi onto the SD card (place it on the root of the SD not in any folder).
2. Place the SD card in the Capture Box and restart the device by unplugging the power cable and then plugging it back in.
3. Once the device starts you will be automatically asked to install the update with the following question: "An update file is available on the SD Card. Do you wish to install?"
4. Use the middle button (Menu button) to navigate to "Confirm" and after that press the left button (Capture button).

NOTE: During the update do not unplug power or remove the SD card. Wait until the installation is complete and the device will automatically restart afterwards.

After the reboot, Capture Box is ready to be used with updated firmware.

---

\*time is calculated based on Full HD images